QUICK START GUIDE

MAS-100 Atmos[®] microbial compressed gas sampler

English - Deutsch - Français - Italiano - Español - РУССКИЙ

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

ENGLISH

Dear customer,

Congratulations on your new MAS-100 Atmos® (Article number: 130.4469/ 11732800001). This quick start guide explains the basics of operating your new instrument - for safe and easy microbial testing of technical and medical gases.

GENERAL INFORMATION

SAFETY SHEET

Carefully read the printout with all safety instructions.

USER MANUAL

To download the full user manual (English only) please visit <https://www.mbv.ch/en/expert-center/downloads/atmos>.

CONTACTS

Manufacturer

MBV AG
www.mbv.ch/service
+41 44 928 30 90
service@mbv.ch

Local service center outside Switzerland

Merck KGaA, Darmstadt, Germany and/or its affiliates
www.sigmaaldrich.com/support/customer-support

TECHNICAL SPECIFICATIONS


Power supply: Input: 100-240V, 1A max, 50-60 Hz (+/- 5%); Output: 5V, 2.1A, 10.5W

Li-ion battery: 10Ah, 7.2V

A complete list of all specifications see user manual chapter 9 «Technical Specifications».

FIRST STEPS

LOCAL USER INTERFACE

CHARGING

- Connect the instrument via USB-C to the power supply
- Check if the integrated LED of the power button is flashing in blue with display switched off or if the flash icon is shown in the upper left corner with display switched on
- Charge until battery is full (continuous blue LED or full battery icon shown)

SWITCH ON/OFF AND STANDBY MODE

- Tap the power button once → Instrument starts up
- Tap the button twice → Instrument switches off
- Tap the power button once while the instrument is switched on → Instrument goes to standby mode and user is logged off

NOTE: While being under pressure the instrument cannot be switched off or go to standby mode. In this case the user is prompted to first decompress the instrument.

FIRST STARTUP: CALIBRATION VALIDITY SEAL

- The instrument comes factory calibrated
- Upon first start of the instrument the calibration interval must be activated by unlocking the calibration validity seal on the instrument screen
- From now on, the calibration validity counts down from 12 months and a warning will be displayed on the screen 1 month before calibration validity ends

CLEANING AND DECONTAMINATION

NOTE: Only perform disinfection when the instrument is at rest. Avoid spray disinfection. Make sure that there is no remaining disinfectant upon starting the instrument.

Housing

Wipe disinfect with one of the following active ingredients:

- Ethanol (EtOH) 70% (v/v)
- 2-Propanol (IPA) 70% (v/v)
- Hydrogen peroxide (H₂O₂) 1%/ acetic acid (AcOH) 5.2%

Perforated lid / dust cover / pressure connection tube / tri-clamp and gasket

- Steam sterilization (autoclaving) at 134°C for 5 mins or 121°C for 20 mins
- Wipe disinfect with 70% Ethanol (v/v) or 70% (v/v) IPA

The list is not exhaustive. For more active ingredients consult the user manual chapter 5 «Care and maintenance».


PERFORM A SAMPLING

The instrument shall only be operated by trained, qualified and authorized personnel. Read all instructions carefully before operating the instrument. Refer to user manual chapter 3 «Operation».

NOTE: It is advised to work with safety goggles and ear protection when handling pressurized gases. Maximum pressure range: 1.2-7 bar absolute (exact range depends on gas type and flow mode see user manual chapter 9 «Technical specifications»). Work in a big and properly ventilated room or use an optional gas exhaust tube.

- Use the pressure connection tube supplied with the instrument to connect the 3/8 inch connector (female) to your gas supply (equipped with a male 3/8 inch connector)
- Fix the 34mm tri-clamp of the pressure connection tube to the sampling head and make sure the connection is tight.
- Open the sampling head (turn counter-clockwise)
- Insert a 90mm nutrient plate aseptically, remove its lid and close the sampling head (turn clockwise)
- Select the sampling parameters of your choice. This pre-selection can be configured (see user manual chapter 4.1.2 «Changing instrument settings».)

- Open the main gas valve (yellow status frame = instrument under pressure)
- Start the sampling. The status of the sampling as well as the current pressure on the instrument are displayed on the screen.
- Close the main valve after the sampling cycle is completed and start the automatic decompression cycle.
- Remove the nutrient plate for incubation and analysis and close the sampling head (turn clockwise)
- Continue with your next sampling location.

ADVANCED FUNCTIONS

ACCESS THE BROWSER-BASED USER INTERFACE (SEE USER MANUAL CHAPTER 2.10.)

- Connect the instrument (USB-C) to a computer (USB-A)
- Open a browser window and enter the IP address **http://192.168.254.1** (no internet needed)
- For login, present the black system administrator hardware key to the key reader and define a personal PIN

INSTRUMENT CONFIGURATION

Refer to chapter 4.1 «Setting up the instrument configuration» in the user manual. Open the submenu 'Instrument settings' and set your preferences.

DATA REPORTING

All data generated by the instrument and by user events are stored in the audit trail of the instrument. Depending on the intended use, the data can be directly transmitted to LIMS or EM systems, viewed in the browser or exported. See all available options in the user manual chapter 3.4 «Data reporting».

OPTIONAL USE CASES: 21 CFR PART 11 SUPPORTING MODE AND DATA INTEGRITY

The instrument offers two optional instrument modes for a use case with or without LIMS or EM system, that can support a workflow that needs to be compliant with 21 CFR Part 11 and needs to meet data integrity requirements (for configuration see user manual chapter 4.3 «21 CFR part 11 supporting instrument modes»).

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The vibrant M, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.