


NOTHING IS EASIER.






HEALTHY INDOOR CLIMATE THANKS TO CLEAN AIR

The MAS-100 NT® 6022 is the certified air sampler for hygiene inspections of HVAC systems. Together with the MAS-100 Eco® it is used to measure microorganisms in the outdoor and supply air.

The aim is to determine the air quality in working and living rooms. And this is done precisely following the latest guidelines according to the hygiene requirements for heating and ventilation (HVAC) systems (SWKI VA104-01/ VDI 6022). That matters, because polluted air endangers health -and reduces physical and mental performance.

MAS-100 NT 6022 and MAS-100 Eco are ergonomic, efficient, and extremely precise. Matching accessories enable maximum measuring reliability and accelerate the measuring process. Because nothing is easier: The MAS-100 hygiene inspection set.







- Rotatable head for floor and lateral ventilation outlets
- Programmable start delay allows measurement on a high tripod
- Extremely robust, durable housing made of anodised aluminium
- Smooth and easy to clean surfaces
- Compatible with all accessories for hygiene inspections



EITHER WAY THE RIGHT DECISION.

MAS-100 ECO

+ YOUR BENEFITS

- Programmable start delay allows measurement on a high tripod
- Extremely robust, durable housing made of anodised aluminium
- Smooth and easy to clean surfaces
- Compatible with all accessories for hygiene inspections

THE IDEAL ACCESSORIES: BECAUSE EVERY DETAIL COUNTS.

TRIPOD

With the tripod the MAS-100 air sampler can be conveyed to the air outlet in no time. It is also compact and lightweight.

QUICK CHANGE ADAPTER

With the quick change adapter the MAS-100 air sampler can easily be attached to the tripod: Fix the adapter plate to the bottom of the air sampler and then click it into the adapter on the tripod. This guarantees a stable mounting and an efficient measuring procedure.

TUBE ADAPTER

The tube adapter made of anodised aluminum fits on all MAS-100 perforated lids and allows the collection of airborne microorganisms at difficult-to-access sampling points such as floor outlets or air ducts.

NBV NEW

AIR FUNNEL

The air funnel supports you in performing induction-free supply air measurements with MAS-100 devices. The two-part system is available with a diameter of 20 or 40cm and can be easily mounted on the air sampler.

HYGIENE INSPECTION CASE

The hygiene inspection case with two rollers can be easily moved from measuring point to measuring point. It protects the device from damage and pollution. At the same time, it offers space for all accessories.

BUY

Invest in a MAS-100 NT 6022 and thus in an extremely robust, durable device that complies with all applicable directives.

DISCOVER

We would be pleased to show you our hygiene inspection set, train you in the correct application, and give you practical tips.

WWW.MBV.CH/HYGIENEINSPECTION

I WOULD BE HAPPY TO ADVISE YOU PERSONALLY.

Corina Keller
Product Manager
corina.keller@mbv.ch
T +41 44 928 30 56



CONTACT US. WE WILL GLADLY HELP YOU WITH A SUITABLE SOLUTION.

Corina Keller
Product Manager
corina.keller@mbv.ch
T +41 44 928 30 56

Dr. Roland Durner
Head of Sales & Marketing
roland.durner@mbv.ch
T +41 44 928 30 82

AIR. NOTHING ELSE.

For more than 30 years, MBV Ltd. has been developing high-quality air samplers to determine the microbial contamination of air. The instruments are used where microorganisms have a negative impact on indoor air quality - and thus harm people or products. In addition to ventilation and air-conditioning systems, for example, in the pharmaceutical, cosmetics, or food industries. Diligence and curiosity have taken us a long way. At least to the top of the world. Learn more at www.mbv.ch.

