

PSD-10 PETRI DISH SAFETY DISPENSER



PSD 10®

Petri dish dispenser

While handling Petri dishes and contact plates is part of the daily routine and does not represent a challenge for trained staff, they often need to be transported to the sampling location and from there to the incubator. The PSD-10 container has been specially developed for the storage and transport of dishes and RODAC plates (using a special adapter). By preventing the dish lid from falling off, it protects the environment from dangerous germs and the dishes from secondary contamination. The Petri dishes remain inside the PSD-10 container during incubation, which further increases protection when handling the dishes inside the incubator.



BENEFITS

- Stackable
- Interchangeable clips in several colours for various growth media or sampling locations
- Absolute protection of valuable samples during transport
- For Petri dishes and RODAC plates
- Autoclavable (as a stack of at least three segments)



SPECIFICATIONS

- Material: polypropylene
- Temperature resistance: 121 °C for 15 minutes (slight shrinkage during the first heat treatment)
- Length 10.3 cm, width 10 cm, height 2.1 cm, weight approx. 38 g (one PSD 10 segment)